# CoviDx™ SARS-CoV-2 Rapid Antigen Test

For in vitro diagnostic use.

Please read this package insert carefully prior to use and strictly follow the instructions.

## INTENDED USE

CoviDx™ SARS-CoV-2 Rapid Antigen Test is a lateral flow immunochromatographic assay intended for the qualitative detection of the nucleocapsid protein antigen from SARS-CoV-2 in nasal (NS) swabs and nasopharyngeal (NP) swabs from patients suspected of a Coronavirus-19 (COVID-19) infection by their health care provider within the first 5 days of symptoms, or from individuals without symptoms or epidemiological reasons to suspect COVID-19 when tested twice over two or three days with at least 24 hours and no more than 36 hours between tests.

CoviDx™ SARS-CoV-2 Rapid Antigen Test aids in the clinical diagnosis of patients suspected of COVID-19. The test is intended for professional use only and should be used in conjunction with other clinical evidence. CoviDx™ SARS-CoV-2 Rapid Antigen Test is for use at the point-of-care and can be used by sites authorized to perform point-of-care testing in near-patient settings.

CoviDx™ SARS-CoV-2 Rapid Antigen Test is not intended for home testing, self-testing or specimen self-collection.

Results are for the identification of SARS-CoV-2 nucleocapsid protein antigen. Antigen is generally detectable in NS and NP samples during the acute phase of infection. Positive results indicate the presence of viral antigens, but clinical correlation with patient history and other diagnostic information is necessary to determine infection status. Positive results do not rule out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease. Positive results are required to be reported to the public health system.

Negative results should be considered presumptive and do not preclude infection. Negative results should not be used as the sole basis for diagnosis, treatment, or other clinical and patient management decisions, including infection control decisions. Negative results should be considered in the context of a patient's recent exposures, history, and the presence of clinical signs and symptoms consistent with COVID-19, and confirmed with a molecular assay, if necessary, for patient management.

CoviDx™ SARS-CoV-2 Rapid Antigen Test does not differentiate between SARS-CoV and SARS-CoV-2.

The performance of this device has not been assessed in a population vaccinated against COVID-19. The performance of the device has not been assessed on specimens from individuals who have been infected with emerging variants of SARS-CoV-2 of public health concern.

## SUMMARY

Coronaviruses are a family of viruses that can cause a range of respiratory illnesses from the common cold to MERS (Middle East Respiratory Syndrome) and SARS (Severe Acute Respiratory Syndrome). SARS-CoV-2 is a new coronavirus that was identified in December 2019 and can cause mild to severe illness (COVID-19). The clinical manifestations include fever, fatigue and other systemic symptoms, which may be accompanied by dry cough and dyspnea which can rapidly develop into severe pneumonia, respiratory failure, acute respiratory distress syndrome, septic shock, multi-organ failure, and severe acid base metabolic disorder. SARS-CoV-2 transmission has been identified primarily through respiratory droplets (sneezing, talking, coughing) and direct contact. Due to its highly contagious nature and global health crises, SARS-CoV-2 has been designated as a pandemic by the World Health Organization (WHO). The ability to rapidly identify COVID-19 patients to enable effective cohorting and prevent the spread of infection is a critical component of combatting the pandemic. Through rapid detection of SARS-CoV-2, CoviDx™ is able to optimize patient isolation decisions and enable accurate cohorting.

## TEST PRINCIPLE

CoviDx™ SARS-CoV-2 Rapid Antigen test is a lateral flow immunochromatographic assay for the rapid detection of nucleocapsid protein antigen from SARS-CoV-2. Colloidal gold particles labeled with anti-SARS-CoV-2 antibody are fixed on the conjugate pad. Anti-SARS-CoV-2 antibody is immobilized on the "T" test line and Goat Anti-Mouse IgG is bound on the "C" control line of the nitrocellulose membrane. SARS-CoV-2 antigen in the sample will bind to anti-SARS-CoV-2 antibody conjugate in the sample pad forming a complex. This complex migrates through the nitrocellulose membrane by capillary action until it reaches the test line where it is captured by the immobilized anti-SARS-CoV-2 antibody, forming a red test line which confirms presence of SARS-CoV-2 antigen. No test line indicates the absence of SARS-CoV-2 antigen. The test contains an internal control (C line), and a red control line will form to indicate that the test was performed correctly and the test reagents are fully functional.

## MATERIALS PROVIDED

<b>25 single use pouched test devices</b>
<b>25 extraction tubes</b>
<b>25 extraction reagent vials</b>
<b>25 sterile swabs</b>
<b>1 workstation</b>
<b>1 package insert</b>
<b>1 Quick Reference Guide (QR)</b>
<b>1 positive control swab (coated with non-infectious recombinant SARS-CoV-2 antigen)</b>
<b>1 untreated negative control swab</b>

## MATERIALS NOT PROVIDED

<b>Timer</b>
<b>Gloves</b>

## KIT STORAGE AND STABILITY

CoviDx™ should be stored at 2-30°C (35.6-86°F) out of direct sunlight. If stored at 2-8°C (35.6-46.4°F), ensure that the test device is brought to room temperature (15-30°C/59-86°F) before testing. Do not freeze. Tests should be used within 1 hour of opening foil pouch. Likewise, extraction reagent vials should be used within one hour of opening. The test device is stable through the expiration date printed on the sealed pouch.

## WARNINGS AND PRECAUTIONS

1. For in vitro diagnostic use only.
2. For prescription use only.
3. Do not use the CoviDx™ test past the expiration date.
4. Use standard precautions for collecting, handling and disposing of patient's nasal or nasopharyngeal swab sample.
5. All samples should be considered potentially hazardous and handled in the same manner as an infectious agent.
6. Use the swab provided in the kit for nasal and nasopharyngeal swab collections.
7. The extraction reagent contains a buffered salt solution (saline). If the extraction reagent comes into contact with the eye or skin immediately flush with copious amounts of water.
8. Adhere to the test instruction to obtain accurate results.
9. Do not remove the test from the foil pouch until ready to use. If the foil pouch is damaged do not use the test.
10. Allow all reagents and samples to reach room temperature before testing 15-30°C (59-86°F).
11. Do not replace the components in this kit with components from other kits.
12. Do not re-use any contents in the kit as they are single-use only.
13. Dispose of used contents as biohazardous wastes in accordance with federal, state, and local requirements.
14. The swab samples should be tested as soon as possible after collection and within 1 hour.
15. The CoviDx™ positive control swab has been prepared from recombinant viral proteins and does not contain infectious material.
16. Do not use the test to evaluate patient samples if either the positive or negative control swab fail to give the expected results.
17. The CoviDx™ test is only intended for nasal and nasopharyngeal swabs samples that are tested directly (i.e. have NOT been placed in transport media).
18. The CoviDx™ test performance has not been validated at altitudes outside of 0-2000 m.

## SAMPLE COLLECTION

CoviDx™ SARS-CoV-2 Rapid Antigen test can be performed on nasal and nasopharyngeal swabs. Use standard procedures for collecting a nasal or nasopharyngeal swab.

### Nasal swab sample collection:

Gently hold and tilt the patient's head with one hand. Carefully insert the swab into the left nostril and gently rotate the swab, ensuring that the swab tip is inserted about 2 cm (< 1 inch) into the nostril until resistance is met at the turbinates. Slowly rotate the swab in a circular path against the nasal wall. Be sure to collect any nasal drainage that may be present on the swab. Gently remove the swab and carefully insert the same swab into the right nostril repeating the collection procedure.

## Nasopharyngeal swab sample collection:

Gently hold the patient's head with one hand, carefully insert the swab into the nostril and slowly enter deep along the bottom of the lower nasal passage. When the top of the swab reaches the back wall of the nasopharyngeal cavity, gently rotate and then slowly remove the swab.


## SAMPLE STORAGE

Freshly collected swab samples should be tested as soon as possible after collection and within 1 hour of collection.

## TEST PROCEDURE

Gloves and suitable eye/face protection should be worn when performing this test. Before testing, read the instructions thoroughly. Allow all reagents and samples to reach room temperature (15-30°C/59-86°F) before testing. Remove the test from the foil pouch just before testing and place on a flat surface.

1. Place an extraction tube in the cardboard workstation. Tear off top of extraction reagent vial and empty the entire contents into the extraction tube.
2. Add the swab to the extraction tube. Rotate the swab approximately 10 times, pressing the swab head against the extraction tube. Squeeze the swab head against the extraction tube to retain as much of the liquid as possible, then remove and discard swab. Immediately proceed to step 3.
3. Place cap on extraction tube, invert the extraction tube, and add 4 drops vertically to the sample well of the test device.
4. Read the test results at 20 minutes. Positive results can be read as early as 15 minutes. Do not read after 30 minutes.

## INTERPRETATION OF RESULTS

### Positive Result

A red line appears in both the test line region 'T' and the control line region 'C'.

A Positive result indicates the presence of SARS-CoV-2 antigen in the sample.


Even if the test line is weak or incomplete it should be interpreted as a positive result. Weak lines may be caused if the antigen concentration is low.

### Negative Result

A red line appears in the control line region 'C' only. No line appears in the test line region "T".

A Negative result indicates the absence of SARS CoV-2 antigen in the sample, or the concentration is below the limit of detection of the test.

### Invalid Result

A control line fails to appear. Insufficient sample volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test device. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

## QUALITY CONTROL


**Internal Control:** CoviDx<sup>®</sup> contains an internal control (C line). A red control line will form to indicate that the test was performed correctly and the test reagents are fully functional. The absence of a red control line indicates an invalid result. Review the procedure and repeat the test with a new test device. If the problem persists, discontinue using the test kit immediately and contact your local distributor or Lumos Diagnostics Technical Support.

**External Control:** External controls may be used to demonstrate that the reagents and assay perform properly. CoviDx™ external controls consist of one (1) positive control swab and one (1) negative control swab. External controls are recommended to be performed with every new shipment, every new lot number, and with every first-time test operator. If the CoviDx™ external control swab does not perform as expected, please repeat the test. If external controls fail on repeat testing, do not perform patient testing and contact your distributor or Lumos Diagnostics Technical Support.

### EXTERNAL CONTROL TEST PROCEDURE

Process each positive and negative control swab exactly the same as a patient swab sample.

1. Place an extraction tube in the cardboard workstation. Tear off top of extraction reagent vial and empty the entire contents into the extraction tube.
2. Add either the positive control swab or the negative control swab to the extraction tube. Rotate the swab approximately 10 times, pressing the swab head against the extraction tube. Squeeze the swab head against the extraction tube to retain as much of the liquid as possible, then remove and discard swab. Immediately proceed to step 3.
3. Place cap on extraction tube, invert the extraction tube, and add 4 drops vertically to the sample well of the test device.
4. Read the test results at 20 minutes. Positive results can be read as early as 15 minutes. Do not read after 30 minutes.

### CLINICAL PERFORMANCE

The performance of the CoviDx™ SARS-CoV-2 Rapid Antigen test was determined by assessing nasal and nasopharyngeal swabs from symptomatic individuals within 5 days of symptom onset who were suspected of COVID-19. The study cohort comprised of males and females, who were racially and ethnically diverse, with a median age of 43 (range 18-83). Of the 227 swabs assessed, 151 were nasal swabs prospectively collected from patients sequentially enrolled by two U.S. outpatient clinics and were representative of a diverse patient population. CoviDx™ SARS-CoV-2 Rapid Antigen testing was performed by nine (9) untrained healthcare professionals who were representative of the intended users and who performed testing in prospectively enrolled patients solely guided by the Quick Reference Instruction (QRI) as well as 1 trained operator who performed testing on retrospectively collected samples. The remaining 76 samples comprised of archived nasal and nasopharyngeal swabs including individuals with variant strains of SARS-CoV-2. The samples were collected from diverse geographic locations, were obtained during high and low prevalence states and were inclusive of newly emerging variant strains including the UK (B.1.1.7.), South African and Denmark variants.

Nasal swabs and NP swabs were collected and handled as described in the CoviDx™ package insert. Prospective samples were tested immediately using the CoviDx™ SARS-CoV-2 Rapid Antigen test. Archived samples were collected and tested on the CoviDx™ SARS-CoV-2 Rapid Antigen test in a randomized and blinded fashion. The performance of the CoviDx™ SARS-CoV-2 Rapid Antigen test was compared to an Emergency Use Authorized molecular (RT-PCR) test for detection of SARS-CoV-2.

Results are presented as the entire cohort as well as stratified by the comparator method cycle threshold (Ct) to estimate the correlation of assay performance in the context of the viral load present in the clinical sample. As presented in the table below, the positive agreement of CoviDx™ SARS-CoV-2 Rapid Antigen test is higher with samples of a Ct count < 35.

The clinical results are summarized in the following tables:

**PROSPECTIVE DATA ONLY (Nasal Swabs Only) - Stratified by CT Value**

	PCR test								
	CT ≤ 30			CT < 35			All (no CT restriction)		
	Positive	Negative	Total	Positive	Negative	Total	Positive	Negative	Total
<b>CoviDx™ SARS-CoV-2 Rapid Antigen Test</b>									
Positive	31	5	36	40	5	45	41	5	46
Negative	0	96	96	0	96	96	9	96	105
Total	31	101	132	40	101	141	50	101	151
Positive Percent Agreement (PPA)	100.0% (95% CI: 89.0% - 100.0%)			100.0% (95% CI: 91.2% - 100.0%)			82.0% (95% CI: 69.2% - 90.2%)		
Negative Percent Agreement (NPA)	95.0% (95% CI: 88.9% - 97.9%)			95.0% (95% CI: 88.9% - 97.9%)			95.0% (95% CI: 88.9% - 97.9%)		

**RETROSPECTIVE DATA ONLY (Nasopharyngeal Swabs Only) - Stratified by CT Value**

	PCR test								
	CT ≤ 30			CT < 35			All (no CT restriction)		
	Positive	Negative	Total	Positive	Negative	Total	Positive	Negative	Total
<b>CoviDx™ SARS-CoV-2 Rapid Antigen Test</b>									
Positive	33	0	33	33	0	33	33	0	33
Negative	0	29	29	3	29	32	4	29	33
Total	33	29	62	36	29	65	37	29	66
Positive Percent Agreement (PPA)	100.0% (95% CI: 89.6% - 100.0%)			91.7% (95% CI: 78.2% - 97.1%)			89.2% (95% CI: 75.3% - 95.7%)		
Negative Percent Agreement (NPA)	100.0% (95% CI: 83.3% - 100.0%)			100.0% (95% CI: 83.3% - 100.0%)			100.0% (95% CI: 83.3% - 100.0%)		

The performance of CoviDx™ has not yet been clinically validated for use in patients without signs and symptoms of respiratory infection or for serial screening applications, and the performance may differ in these populations.

Clinical studies in asymptomatic patients using serial testing are ongoing to establish clinical performance.

## ANALYTICAL PERFORMANCE

### Microbial Organism Cross Reactivity/Interference

No cross-reactivity or interference was observed for the following microorganisms at the listed concentrations.

Potential Cross-Reactant/Interferent	Test Concentration CFU/mL or PFU/mL	Reviewer Result
Adenovirus (Type 7A Species B)	1×10 <sup>5</sup> PFU/mL	Negative
Human Metapneumovirus (hMPV) 16 Type A1	1×10 <sup>5</sup> PFU/mL	Negative
Parainfluenza virus 1-4 (Type 1)	1×10 <sup>5</sup> PFU/mL	Negative
Parainfluenza virus 1-4 (Type 2)	1×10 <sup>5</sup> PFU/mL	Negative
Parainfluenza virus 1-4 (Type 3)	1×10 <sup>5</sup> PFU/mL	Negative
Parainfluenza virus 1-4 (Type 4 B)	1×10 <sup>5</sup> PFU/mL	Negative
Influenza A H1N1	4.9×10 <sup>4</sup> PFU/mL	Negative
Influenza B Virus	4.1×10 <sup>4</sup> PFU/mL	Negative
Enterovirus (Type 68)	1×10 <sup>5</sup> PFU/mL	Negative
Respiratory syncytial virus (Type A)	1×10 <sup>5</sup> PFU/mL	Negative
Rhinovirus	4.1×10 <sup>4</sup> PFU/mL	Negative
<i>Haemophilus influenzae</i> (Type B) Eagan	1×10 <sup>6</sup> CFU/mL	Negative
<i>Streptococcus pneumoniae</i>	1×10 <sup>6</sup> CFU/mL	Negative
<i>Streptococcus pyogenes</i>	1×10 <sup>6</sup> CFU/mL	Negative
<i>Candida albicans</i>	1×10 <sup>6</sup> CFU/mL	Negative
Pooled human nasal wash – representative of normal respiratory microbial flora	N/A	Negative
<i>Bordetella pertussis</i>	1×10 <sup>6</sup> CFU/mL	Negative
<i>Mycoplasma pneumoniae</i>	1×10 <sup>6</sup> CFU/mL	Negative
<i>Chlamydia pneumoniae</i>	1×10 <sup>6</sup> CFU/mL	Negative
<i>Legionella pneumophila</i>	1×10 <sup>6</sup> CFU/mL	Negative
<i>Staphylococcus aureus</i>	1×10 <sup>6</sup> CFU/mL	Negative
<i>Staphylococcus epidermidis</i>	1×10 <sup>6</sup> CFU/mL	Negative
SARS-coronavirus	Positive Detection by PCR	Negative
Human coronavirus 229E	2.0×10 <sup>4</sup> PFU/mL	Negative
Human coronavirus OC43	1×10 <sup>5</sup> PFU/mL	Negative
Human coronavirus NL63	4.9×10 <sup>4</sup> PFU/mL	Negative
MERS-coronavirus	4.1×10 <sup>4</sup> PFU/mL	Negative

## INTERFERING SUBSTANCES

None of the tested substances exhibited any interference to CovidX™ within the given concentrations.

Substance	Test #	Reviewer	True Value
Whole Blood (4%)	1	POS	POS
	1	POS	POS
	1	POS	POS
Naso GEL (NeilMed) (5% v/v)	2	POS	POS
	2	POS	POS
	2	POS	POS
CVS Nasal Drops (Phenylephrine) (15% v/v)	3	POS	POS
	3	POS	POS
	3	POS	POS
Afrin (Oxymetazoline) (15% v/v)	4	POS	POS
	4	POS	POS
	4	POS	POS
CVS Nasal Spray (Cromolyn) (15% v/v)	5	POS	POS
	5	POS	POS
	5	POS	POS
Zicam (5% v/v)	6	POS	POS
	6	POS	POS
	6	POS	POS
Homeopathic (Alkalol) (1:10 dilution)	7	POS	POS
	7	POS	POS
	7	POS	POS
Sore Throat Phenol Spray (15% v/v)	8	POS	POS
	8	POS	POS
	8	POS	POS
Fluticasone Propionate (5% v/v)	9	POS	POS
	9	POS	POS
	9	POS	POS
Chloraseptic (Menthol/Benzocaine) (1.5 mg/mL)	10	POS	POS
	10	POS	POS
	10	POS	POS
Tobramycin (4 µg/mL)	11	POS	POS
	11	POS	POS
	11	POS	POS
Mucin (0.5%)	12	POS	POS
	12	POS	POS
	12	POS	POS
Mupirocin (10 mg/mL)	13	POS	POS
	13	POS	POS
	13	POS	POS
Tamiflu (Oseltamivir Phosphate) (5 mg/mL)	14	POS	POS
	14	POS	POS
	14	POS	POS

## HOOK EFFECT

No high dose hook effect was observed up to  $9 \times 10^5$  TCID<sub>50</sub>/mL of heat-inactivated SARS-CoV-2 with CovidX™.

## LIMIT OF DETECTION

The limit of detection of CovidX™ was determined to be  $9.5 \times 10^2$  TCID<sub>50</sub>/mL.

## LIMITATIONS OF THE TEST

1. CoviDx™ can only be used for the qualitative detection of SARS-CoV-2 antigen and cannot determine the quantity of SARS-CoV-2 antigens in samples.
2. CoviDx™ should be used for the detection of SARS-CoV-2 antigen in nasal swab and nasopharyngeal swab samples only.
3. Results from antigen testing should not be used as the sole basis to diagnose or exclude SARS-CoV-2 infection. All results must be interpreted together with other clinical information including symptoms/signs, medical history, and other laboratory tests.
4. If the CoviDx™ test result is negative and clinical symptoms persist, it is recommended to repeat sampling or use other testing methods for testing. A negative result cannot preclude the possibility of exposure to or infection with SARS-CoV-2 virus at any time.
5. Inadequate sample collection, transportation and processing, or low virus titer in the sample, may lead to false negative results.
6. The performance of CoviDx™ has not yet been evaluated for use in patients without signs and symptoms of respiratory infection or for serial screening applications, and the performance may differ in these populations.
7. Clinical studies in asymptomatic patients using serial testing are ongoing to establish clinical performance.

## TECHNICAL SUPPORT

Please contact Lumos Diagnostics Technical Support at 941.556.1850 for assistance.

## AVAILABILITY

CoviDx™ Rapid Antigen Test Part# CP0033

## REFERENCES

1. Weiss SR, Leibowitz JL. Coronavirus pathogenesis. *Adv Virus Res.* 2011; 81:85-164.
2. Masters PS, Perlman S. Coronaviridae. In: Knipe DM, Howley PM, eds. *Fields Virology*. 6th ed. Philadelphia, PA: Lippincott Williams & Wilkins; 2013:825-58.
3. Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses. *Trends Microbiol.* 2016; 24:490-502.
4. Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses. *Nat Rev Microbiol.* 2019; 17:181-192.
5. Wong G, Liu W, Liu Y, Zhou B, Bi Y, Gao GF. MERS, SARS, and Ebola: the role of super-spreaders in infectious disease. *Cell Host Microbe.* 2015; 18:398-401.
6. <https://covid19.who.int/>

 <b>Manufacturer</b>	 <b>Date of manufacture</b>	 <b>In vitro diagnostics</b>
 <b>Caution</b>	 <b>Consult instructions for use</b>	 <b>Catalogue number</b>
 <b>Temperature limit</b>	 <b>Number of tests</b>	 <b>Batch code</b>
 <b>Do not re-use</b>	 <b>Sterilized by ethylene oxide</b>	 <b>Negative control</b>
 <b>Expiration date</b>	 <b>Do not use if package is opened or damaged</b>	 <b>Positive control</b>

### Manufacturer

 Rapid Pathogen Screening, Inc.

7040 Professional Parkway, Suite B, Sarasota, FL 34240 USA

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